Computer Vision Techniques

Neuralink, US 2020/0085508 A1, 3/19/2020 4/28/2024

Background

- Conventional surgical techniques prohibitively expensive, limit BCI innovation (cannot sequentially implant large # of electrodes into brain)
- Conventional robotic surgery inferior to human surgeons in terms of precisely positioning implants, target tissues, insertion needles (due to glare, lighting, reflective elements)
- Robotic surgery system use fluorescence of certain elements, special lighting, computer vision techniques for bio-compatible electrodes in biological tissue using robotic assemblies


Outline

• Goal: Discuss computer vision techniques in connection with robotic surgery system (implantable device engagement + targeting + insertion verification subsystems)

Outline

- Embodiments:
 - Target Components:
 - System for implanting a device Front & Side Views
 - Inserter head with needle for inserting probe
 - Implantable Device Engagement Components:
 - Probe device stage with cartridge-pillbox assembly
 - Needle and engagement component illuminated by white light
 - Engagement component fluorescing in response to irradiation
 - Needle and Engagement Component taken using Red Light
 - Electrode Implantation:
 - Insertion Needle and Pincher on Needle Pincher Cartridge
 - Implantation of Electrodes in a Target Tissue Proxy Substance
 - Example of Electrodes Implanted in Brain Tissue
 - Insertion Verification Components:
 - Example Verification Components of a System for Robotic Surgical Implantation
 - Robotic Surgical Implantation Techniques:
 - Flowcharts for Device Engagement and Robotic Surgical Implantation
 - Robotic Surgical Implantation (including Targeting, Verification)
 - Example Computing System for Robotic Surgery guided by Computer Vision
 - Example Components of Computing System for Robotic Surgery guided by Computer Vision

Robotic System for implanting a device – Front View



Robotic System for implanting a device - Side View






FIG. 1B

Inserter head with needle for inserting probe



Probe device stage with cartridge-pillbox assembly



Needle and Engagement Component illuminated by white light

- Both needle & engagement component appear blurry
- Difficult for system to discern loop of engagement component
- Needle appears dimly lit
- May not be able to determine positions/orientations of needle & engagement component to engage them
- Prone to erroneous detections of needle & engagement component




FIG. 4

Engagement component fluorescing green light in response to irradiation

- Engagement component fluoresces green light in response to irradiation at near-UV wavelengths (405nm)
- Appears much more sharper than in figure 4 system can locate engagement component more accurately
- Color data may help identifying the engagement component




FIG. 5

Needle and Engagement Component taken using Red Light

- Red light provdes technical advantage reflects more strongly and clearly from metal with needle
- Appears **sharper** than in figure 4
- Red light with red backing enable cameras to resolve images with very defined edges on needle tip
- Uses red light to detect position/orientation of needle
- Avoids erroneous detections, false positives, and false negatives
- Note: images can be in color (may be helpful for identification)





FIG. 6

Insertion Needle and Pincher on Needle Pincher Cartridge (NPC)


• NPC:

- Portion of inserter head that makes direct contact with brain
- Can be replaced mid-surgery in under minute
- Milled from $40 \mu m$ diameter tungstenrhenium wire-stock electrochemically etched to 24m diameter along inserted length



Implantation Process of Electrodes in a Target Tissue Proxy Substance

- Needle inserts first thread (holds plurality of electrodes)
- Insert second thread with second plurality of electrodes



Example of Electrodes Implanted in Brain Tissue

- 96 polymer threads, 32 electrodes/thread
 - 5-50 μm thread width, **nominal thread** thickness 4-6 μm , length ~20mm
- Thin Hair: $60\mu m$ thick (12-15x thicker than threads)
- Small size, greater flexibility, greater biocompatibility
- Probes can remain implanted for long periods of time without triggering immune responses, minimize tissue displacement





FIG. 9

Example Verification Components of a System for Robotic Surgical Implantation




- process images from cameras via CV heuristic to determine implantable electrode device and/or insertion needle
- send further instructions to robotic assemblies (ex. correct positioning implantation of electrode device

no further motions needed, robotic surgery proceed to next stage)


Device Engagement and Robotic Surgical Implantation

Note: The following diagrams may involve more components (e.g., cameras, illumination components, computing system, etc.) but are not shown for simplicity





Step 1:

IRRADIATE POLYMER PORTION OF IMPLANTABLE DEVICE WITH NEAR UV LIGHT $\underline{1102}$




Step 2:






Step 3:


IMAGE POLYMER PORTION WITH SECOND CAMERA ALONG SECOND DIRECTION 1106






Step 4:

PROCESS IMAGES TO TRIANGULATE 3D LOCATIONS AND/OR ORIENTATIONS OF IMPLANTABLE DEVICE AND INSERTION NEEDLE IN SURGICAL FIELD 1108




Step 5:

ILLUMINATE AN INSERTION NEEDLE USING VISIBLE LIGHT




Step 6:

IMAGE INSERTION NEEDLE WITH FIRST CAMERA 1112




Note: real image not available for show

Step 7: ROBOTICALLY ENGAGE THE POLYMER PORTION OF THE ELECTRODE WITH THE INSERTION NEEDLE 1114 needle assembly




IMPLANT INSERTION NEEDLE AND ELECTRODE 1118


Step 9:




Robotic Surgical Implantation



Step 1:




POSITION PILLBOX CARTRIDGE ASSEMBLY AND NEEDLE IN SURGICAL FIELD Step 2: 1204 pillbox needle cartridge assembly assembly



Step 4:

OPTIONALLY I MAGE TARGET SURGICAL TISSUE 1208






Note: real image not available for show


Step 5:


contour location (processed with vectors, contours, etc.)



Step 6:

IMPLANT IMPLANTABLE DEVICE AND/OR NEEDLE INTO BIOLOGICAL TISSUE 1212





Step 8:




FIG. 10

Verification: computing system determines whether to end system or return to step 4 for further correction





FIG. 9


Targeting during Robotic Surgical

Implantation APPLY LIGHT TO SURGICAL FIELD IMAGE SURGICAL FIELD TRANSMIT IMAGE TO COMPUTING SYSTEM DETERMINE ONE OR MORE TARGET TISSUE SITES FOR IMPLANTATION 1202 PROCESS IMAGE TO DISTINGUISH BIOLOGICAL STRUCTURES, TISSUE DISPLAY PROCESSED IMAGE 1310 GENERATE IMPLANTATION TARGET BASED ON PROCESSED IMAGE GENERATE SURGICAL PLAN BASED ON IMPLANTATION TARGETS. **OBTAINED IMAGE** FIG. 13

Step 1:

APPLY LIGHT TO SURGICAL FIELD 1302






Step 3:


TRANSMIT IMAGE TO COMPUTING SYSTEM $\underline{1306}$






Step 4:

PROCESS IMAGE TO DISTINGUISH BIOLOGICAL STRUCTURES, TISSUE 1308




Step 5:

DISPLAY PROCESSED IMAGE 1310



Step 6:


GENERATE IMPLANTATION TARGET BASED ON PROCESSED IMAGE 1312





GENERATE SURGICAL PLAN BASED ON IMPLANTATION TARGETS, **OBTAINED IMAGE** Step 7: 1314

(in computer language, of course)

Verification during Robotic Surgical Implantation







Step 2:

TRANSMIT IMAGE TO COMPUTING SYSTEM 1404






FIG. 9




Step 3:

PROCESS OBTAINED IMAGE ACCORDING TO COMPUTER VISION HEURISTIC TO IDENTIFY ELECTRODE DEVICE, INSERTION NEEDLE, TARGET SITE 1406




DETERMINE WHETHER ELECTRODE DEVICE IS PROPERLY IMPLANTED $\underline{1408}$


Step 4:




Example Computing System for Robotic Surgery guided by Computer Vision



Example Components of Computing System for Robotic Surgery guided by Computer Vision Vision




Video: Implanting Thread Process



Source: https://zhuanlan.zhihu.com/p/574574159

Video: Implanting Thread Process



Source: https://zhuanlan.zhihu.com/p/574574159